

NEW

Hospital TV Remote

CLEAN REMOTE®

Clean Remote® is scientifically tested to carry 99% Less Bacteria than other remotes

In several studies, TV remotes have been revealed to be a leading carrier of bacteria in a patient's room. Not only are they handled continuously throughout the day by patients, but remotes are also handled by staff and visitors. Standard remotes and nurse call pillow remotes contain many cracks and crevices, which make it impossible to properly disinfect. The easy-to-wipe non-porous surface of the Clean Remote® makes it simple to clean and disinfect.

- **Simple One Touch Set-up** (no codes needed)
- **Works Over 200 Brands of Commercial and Retail TVs**
- **No Menu Button** (patients and guests cannot change settings)
- **Hidden Battery Compartment** (helps prevent batteries from disappearing)
- **Uses 2 AAA batteries**

The Clean Remote® offers direct channel access (numbers 1 - 9 and 0) whereas most pillow remotes have only channel up and channel down arrows, which patients cycle through all the channels to get back to a specific channel.

Hidden battery compartment and patent pending hinged battery door helps eliminate lost battery doors and deters patients from removing batteries.



Note: Clean Remote® will not operate with Lodgenet or On Command movie systems. Also will not work TV/DVD, VCR or DVD/VCR combination units.



Model 1234



Call Today For More Information Or To Order

Toll Free 1-800-288-7351

Toll Free Fax 1-888-691-5268